

## TF513A/C/D Circular Sample Cutter

### Application

Fabric Circular Sample Cutter / GSM Cutter. The circular sample cutter cuts fabric sample in various dimension, such as 100 cm<sup>2</sup> for GSM, Dia. 140 mm for Martindale test, etc. Cut accurate fabric circular samples with smooth edges by drawing action even difficult materials such as fabric, thin films, tissue paper, and corrugated cardboard and synthetic leather can be cut by the sample cutter conveniently.

The cutting pad is made of porous rubber, which allows the sample cutter to cut into the base with complete safety.

### Models

- TF513A Circular Sample Cutter -- Dia. 113 mm / 100 cm<sup>2</sup> --For European Fabric Yield System
- TF513C Circular Sample Cutter -- Dia. 38 mm / 12 cm<sup>2</sup> -- For Martindale Abrasion and Pilling Tester
- TF513D Circular Sample Cutter Dia. 140 mm / 154 cm<sup>2</sup> --For Martindale Abrasion and Pilling tester

**Weight** A: 3 kg, C: 1 kg

**Dimensions** A: 175 x 170 x 130 mm (L x W x H)

C: 115 x 110 x 150 mm (L x W x H)



TF513A



TF513C



TF513D

**Signature**

TESTEX INSTRUMENT LTD

Vivian Fu



Our Instruments

[www.Testextextile.com](http://www.Testextextile.com)

Your Success