


Certificate of Conformity

SDC Phenolic Yellowing Test Papers

Product Code: 4225

Unit Size	100mm x 75mm - 50 Piece pack	
Batch:	YP9-081005	
Year of Manufacture:	2018	

Statement of Conformity

This is to certify that Phenolic Yellowing Test Papers have been produced and independently tested to conform to the specifications of BS EN ISO 105 – X18:2007. M&S C20B and NEXT TM43

Independent Testing

Where applicable, to ensure that all results are accurate, unbiased and impartial, testing is carried out in **Independent UKAS Accredited Test Laboratories** and sampling is conducted using ISO guidelines as a minimum requirement.

In addition to performance testing multiple samples of each batch of impregnated paper is subjected to chemical analysis to determine the level of active ingredient in the paper. These results are checked versus internal specifications and the level set out in clause 7.1 of BS EN ISO 105 – X18:2007.

Packaging

Please ensure that the product is retained in its original SDCE packaging, which contains the correct batch details. All product packaging is tested to ensure there are no contaminants or BHT content that could cause fabric yellowing. Pack sizes are measured and approved using UK weights and measures act.

Signed:	
	J Foister – Technical Manager