

TY360 Wrap Reel

Wrap Reel, is to produce skeins of yarn of a pre-determined length and number of turns for count and strength testing. 1 Meter, 36" or 54" circumference collapsible swift (specify). Wrap reel is completed with yarn package stand and pre-tension device, and fitted with pre-determined counter.

Specifications

- Circumference of winch 1000 ±1 mm
- Number of wraps 1 ~ 9,999 adjustable
- Pre-tension 2 ~ 100cn
- Traveling reciprocating distance 35 mm
- Spacing of spindles 60 mm
- Reel speed 20 ~ 280 r/min (variable)

Standards

ISO 2060, ASTM D1907/2260, BS 2010, DIN 53830, GB/T4743, GB/T14343, GB/T6838



Power 220 / 110 V 50 / 60 Hz

Weight 75 kg

Dimensions 960 x 760 x 670 mm (L x W x H)

Signature

TESTEX INSTRUMENT LTD



Our Instruments

www.Testextextile.com

Your Success